

# Experiment 29

## Infra-Red Spectroscopy

Jay C. McLaughlin  
Colorado Northwestern Community College

Name:

CC-BY-SA - March 26, 2022

Date:

### Key Objectives

1. Basic understanding of interaction of infra-red (IR) radiation and chemical bonds.
2. Create a flow chart to identify organic molecules using IR spectroscopy.
3. Identify functional groups in organic molecules.

### Discussion

#### A. Infrared Spectroscopy Background

We recall (hopefully!!) from last semester that electrons around an atom can absorb and emit light of only specific colors (line spectra) due to the quantized nature of the electron orbitals. Each compound tested had a unique spectrum, and we used the intensity of the peak to give us information about the concentration of the compound using Beer's Law.

Figure 29.1 below shows the entire Electromagnetic Spectrum of which in this lab we will obviously focus on the IR region which is located to the right of the Visible region, and therefore is of lower energy, and longer wavelength.




Figure 29.1: The Electromagnetic Spectrum. credit: "The Electromagnetic Spectrum" by LibreTexts is licensed under notset.

Scientists often like to use units which provide nice numbers to work with, thus for IR radiation we generally use the wavenumber ( $\nu$ ) instead of the wavelength ( $\lambda$ ) of the radiation. The relationship between the two is shown in the equation below. This makes the mid IR range between 400-4000

## Experiment 29 IR Spectroscopy

$\text{cm}^{-1}$ . Note that the wavenumber is proportional to energy of the radiation, as wavenumber increases so does the energy of the radiation.

$$\nu = \frac{1}{\lambda} \text{ where } \nu \text{ is measured in units of } \text{cm}^{-1}$$

The wavelength of the electromagnetic radiation determines how it interacts with a molecule. For the visible radiation we saw that the interaction correlated with the transitions of electrons between orbitals and gave rise to Line Spectra (and the color of compounds) that were unique for each element. IR radiation interacts with molecules primarily by being converted into different molecular vibrations, causing chemical bonds to stretch and bend.




Figure 29.2: How Electromagnetic Radiation Interacts with Molecules. credit: unknown

Infrared radiation interacts with the chemical bonds between atoms, and produces a unique spectrum which can tell us information about the type of bonds in a compound. Our normal picture of a molecule (similar to the Lewis Structures/Models build in the previous lab) used rigid bonds between the atoms. This picture is a bit simplified, and instead of rigid bonds between atoms we should picture the bonds as a spring which can vibrate, stretch, and bend in many different ways. Infrared radiation causes vibrations of the atoms which cause the molecular bonds to stretch or bend. There are several types of bending and stretching motions illustrated in Figure 29.3 below (for an dynamic picture see [http://en.wikipedia.org/wiki/Infrared\\_spectroscopy](http://en.wikipedia.org/wiki/Infrared_spectroscopy)). As a molecule becomes larger and more complicated the number of vibrational modes (or degrees of freedom) increases as  $3N-5$  for linear molecules and  $3N-6$  for non-linear molecules, where  $N$  = number of atoms. Obviously very large molecules can give rise to complicated spectra.




Figure 29.3: Basic Bending and Stretching Modes for a Molecule. Animations located at [https://en.wikipedia.org/wiki/Infrared\\_spectroscopy](https://en.wikipedia.org/wiki/Infrared_spectroscopy). credit: unknown

A typical IR spectrum can be seen below in (Figure 29.4 and may seem incredibly complicated at first due to the large number of peaks compared to the visible spectra we took last semester. However, all is not lost, we will focus on learning only a few of the basic rules; enough to identify the functional group in most cases. The very complex looking region below  $1500\text{ cm}^{-1}$  is called the **Fingerprint Region** and can mostly be ignored. The fingerprint region contains bond stretches between backbone atoms (C-C mostly) which almost all organic molecules have, and there are often too many overlapping absorptions to be any help. Instead, concentrate on the **Functional Group Region** between  $1500\text{--}4000\text{ cm}^{-1}$ . This region contains peaks that correspond to stretches in functional groups, and is all that is usually needed to determine what type of compound we have. Table ?? on the next page lists each functional group and the most common absorption frequencies for each. By comparing between Table ?? and each IR spectra the functional group for each molecule can be determined and a model of the molecule can begin to take shape. In order to fully determine the structure of a molecule, it is generally necessary to have data from other techniques, such as chemical formula from a Mass Spectrometer (MS - Lab 32) or number and type of carbons and hydrogens from Nuclear Magnetic Resonance (NMR - Lab 33), to predict or confirm a specific structure.




Figure 29.4: Example IR Spectrum. credit: <https://www.chem.ualgary.ca/courses/351/Carey5th/Ch13/ch13-ir-1.html>

Organic chemists have observed and catalogued the types and locations of Infrared absorptions produced by variety of chemical bonds. The data is generally presented through Tables or Figures linking functional groups with specific regions of the Infrared Spectrum. These are general rules only, as the exact bonding, and nearby functional groups can all influence what frequencies a molecules will absorb. With the advent of computers, much of this information is online, and can even be searched to find a matching spectrum. We, however, will practice matching up spectra and functional groups the old fashioned way in order to gain a basic understanding of the technique.

## Experiment 29 IR Spectroscopy

Graphically the different regions can be seen in Figure 29.5 and 29.6 below. Specific regions and functional groups are discussed in more detail in several sections below.




Figure 29.5: Graphical Representation of Typical IR Frequencies. credit: InfraredSpectraofSomeCommonFunctionalGroups. (2019, September 3). <https://chem.libretexts.org/@go/page/183009>




Figure 29.6: Graphical Representation of Typical IR Frequencies. credit: unknown

### Alkanes, Alkenes, Alkynes and Aromatics

The absorption of the alkane part of a molecule consist primarily of C-H stretching and bending. When other functional groups are present these absorptions are often difficult to distinguish. It is often the lack of absorptions from other functional groups that leads one to identify a molecule as containing only an alkane functional group.

Functional Group	Wavenumber (cm <sup>-1</sup> )	Intensity	Comments
Alkane			
C-H stretch -C-H bend	2850-3000 1350-1480	strong variable	
Alkene			
C=C stretch	1620-1680	variable (weak)	*Only alkenes and Aromatics show C-H absorption > 3000
=C-H stretch	3000-3100	medium	
=C-H bend	650-1000	strong	
Alkyne			
C≡C stretch	2100-2300	variable	*terminal only (can be obscured by other groups (OH)) *terminal only
C≡C-H stretch	3300	strong	
C≡C-H bend	600-700	strong	
Aromatic			
=C-H stretch	3000-3100	medium	*Only alkenes and Aromatics show C-H absorption > 3000
C-C stretch	1400-1600	medium, sharp	multiple bands (1400-1500 and 1585-1600)
overtones	1650-2000	weak	reflect substitution pattern on ring (ortho, meta, para)
C-H "oop"	675-900	weak	out-of-plane C-H bends reflects substitution patterns

Table 29.1: Alkanes, Alkenes, and Alkynes

The IR spectrum of octane is shown in Figure 29.7 below. Note the strong bands in the 3000-2850 cm<sup>-1</sup> region due to C-H stretch. The C-H scissoring (1470 cm<sup>-1</sup>), methyl rock (1383 cm<sup>-1</sup>), and long-chain methyl rock (728 cm<sup>-1</sup>) are noted on this spectrum. Since most organic compounds have these features, these C-H vibrations are usually not noted when interpreting a routine IR spectrum.

The addition of double bonds (C=C) gives rise to a new diagnostic peak, and shifts the frequency of the C=C-H stretch to a high value than the C-C-H stretch, generally slightly greater than 3000 cm<sup>-1</sup>. This is a very useful diagnostic fact since only alkenes and aromatics show a C-H stretch > 3000 cm<sup>-1</sup>.

## Experiment 29 IR Spectroscopy




Figure 29.7: IR spectra of Octane and 1-octene. credit: credit: InfraredSpectraofSomeCommonFunctionalGroups. (2019, September 3). <https://chem.libretexts.org/@go/page/183009>

Alkynes triple bonds ( $C\equiv C$ ) provide additional features that can be identified. Typical bands are shown in Table 29.8. Triple bonds ( $C\equiv C$  and  $C\equiv N$ ) have a distinctive (sharp, though often small) absorption at  $2200\text{ cm}^{-1}$ . Distinguishing between the two usually requires additional information, a terminal  $C\equiv C-H$  bond will have a peak at  $3270\text{--}3330\text{ cm}^{-1}$  due to the C-H stretch, and a peak at  $610\text{--}700\text{ cm}^{-1}$  due to C-H bend.

1-hexyne illustrates the typical features found in a compound with a triple bond (peak at  $2260\text{--}2100\text{ cm}^{-1}$ ), and also illustrates the terminal  $C\equiv C-H$  stretch for when the triple bond is located between 2nd and 3rd carbons.




Figure 29.8: IR spectra of 1-hexyne and Toluene. credit: credit: InfraredSpectraofSomeCommonFunctionalGroups. (2019, September 3). <https://chem.libretexts.org/@go/page/183009>

The spectrum of toluene is shown in Figure 29.8. The  $=C-H$  stretches of aromatics ( $3100$ ,  $3070$ ,  $3030$ ) and the  $C-H$  stretches of the alkyl (methyl) group ( $2925$ ). The characteristic overtones are seen from about  $1650\text{--}2000$ . Also note the carbon-carbon stretches in the aromatic ring ( $1620$ ,  $1510$ ,  $1470$ ), the in-plane  $C-H$  bending ( $1090$ ,  $1040$ ), and the  $C-H$  out-of-plane ( $740$ ).

## Aldehydes, Ketones, Carboxylic Acids and Esters

The typical frequencies for double bonds (C=C, C=O, and C=N) overlaps significantly and it is generally necessary to have additional information to identify them exactly. However, in this class, we can exclude C=N bonds and in conjunction with other information given identify them.

Functional Group	Wavenumber (cm <sup>-1</sup> )	Intensity	Comments
Aldehyde			
C=O	1720-1740	strong	
terminal C-H	2720-2750 and 2820-2850	medium	Two absorptions
Ketone			
C=O	1680-1750	strong/sharp	Lack of other bonds normally indicates a ketone
Carboxylic Acid			
C=O	1710-1780	strong	
O-H stretch	2500-3300	strong, broad	
O-H bend	900-955 and 1400-1450	medium	
C-O stretch	1210-1320	strong	
Ester			
C=O stretch	1735-1750	strong	
C-O stretch	1000-1300		Two or more absorptions
Amides			
C=O stretch	1630-1690	strong	
N-H stretch	3100-3500	medium	
N-H bend	1550-1640	medium	unsubstituted = 2 absorptions

Table 29.2: Double Bonds




Figure 29.9: IR spectra Butanal and 2-butanone. credit: InfraredSpectraofSomeCommonFunctionalGroups. (2019, September 3). <https://chem.libretexts.org/@go/page/183009>



## Experiment 29 IR Spectroscopy




Figure 29.10: IR spectra Hexanoic Acid and Ethyl Ethanoate. credit: <https://orgchemboulder.com/Spectroscopy/irtutor/tutorial.shtml>




Figure 29.11: IR spectra for Butanamide and a Pufferfish. credit: Left - NIST Chemistry WebBook <http://webbook.nist.gov/chemistry> Right - <https://clipartix.com/clipart-fish-image-13780/>

Carbonyl Groups (C=O stretches) typically have a sharp finger around 1715 for a ketone. The lack of additional features is indicative of a ketone.



Carboxylic acids show a strong, wide band for the O–H stretch. Unlike the O–H stretch band observed in alcohols, the carboxylic acid O–H stretch appears as a very broad band in the region 2500-3300, centered at about 3000. This is in the same region as the C–H stretching bands of both alkyl and aromatic groups. Thus a carboxylic acid shows a "messy" absorption pattern in the region 2500-3300, with the broad O–H band superimposed on the sharp C–H stretching bands. The reason that the O-H stretch band of carboxylic acids is so broad is because carboxylic acids usually exist as hydrogen-bonded dimers.



## Alcohols, Ethers, and Amines

Functional Group	Wavenumber (cm <sup>-1</sup> )	Intensity	Comments
Alcohol			
O–H stretch	3200-3500	broad, strong	H-bonding present
O–H stretch	3500-3700	strong, sharp	Free O-H (no H-bonding) (sometimes missing)
C–O stretch	1000-1300		Two or more absorptions
Ether			
C–O	1000-1300	strong	Lack of other absorptions
Amine			
N–H stretch	3000-3300	medium, sharp	Same region as O–H but weaker and sharper. 1° = 2 absorptions, 2° = 1 absorption (weak), and 3° = none
N–H bending	1550-1650	medium	1° only, 2° rarely seen, ° = none, often difficult to tell due to being near C=O region.
N–H wag	650-900	strong, sharp	1° and 2° only
C–N stretch	1050-1350	medium-weak	

Table 29.3: Alcohols, Ethers and Amines

Figure 29.12: IR spectra Ethanol and Ethoxyethane. credit Left - <https://orgchemboulder.com/Spectroscopy/irtutor/tutorial.shtml> right - unknownFigure 29.13: IR spectra for Aniline (1°) and Ethylethanamine(2°). credit: <https://orgchemboulder.com/Spectroscopy/irtutor/tutorial.shtml>

## Experiment 29 IR Spectroscopy




Figure 29.14: IR spectra for N,N-diethylethanamine (3°) and a Platypus. credit: left - <https://orgchemboulder.com/Spectroscopy/irtutor/tutorial.shtml> right - Public Domain

### 3200-3500 cm<sup>-1</sup> Region

This region has several specific functional groups (Alcohols, Amines (1° and 2°) and Carboxylic Acids) which have characteristic peaks in this region. It is the shape of the peak, more than the position which gives clues as to the appropriate functional groups. If the region contains no absorption then these functional groups are likely missing.




Figure 29.15: IR Identification in the 3500-3200 Region. credit - <https://toeholds.wordpress.com/2010/09/21/pictorial-guide-to-interpreting-infrared-spectra/>

## How to Analyze an IR Spectrum

The task of analyzing an IR spectrum can seem daunting as there are many functional groups and many different bending and stretching combinations present. There are many different approaches one may take to analyze a spectrum, the following is just one example.

Note that just as often as a specific absorption indicates a class of compounds, the lack of a specific absorption can also help exclude compounds from consideration.

1. Examine the  $1700 \pm 1000 \text{ cm}^{-1}$  carbonyl region. If a carbonyl group is present, look for the different secondary absorptions associated with each specific functional group to verify your assignment.
  - (a) Aldehyde: C–H sharp absorption at 2720 and 2820  $\text{cm}^{-1}$  (sometimes seen as shoulders)
  - (b) Ketone: no additional absorptions (i.e. lack of Aldehyde/Carboxylic Acid/Ester/Amide peaks)
  - (c) Carboxylic Acid: broad O–H absorption 3000-3500  $\text{cm}^{-1}$  + others
  - (d) Ester: C–O absorption 1000-1300  $\text{cm}^{-1}$  (1 or more)
  - (e) Amide: broad N–H absorption 3000-3500  $\text{cm}^{-1}$  + others
2. Examine the 3000-3500  $\text{cm}^{-1}$  region for broad features (not the C-H stretches).
  - (a) Alcohol: broad, smooth absorption, confirmation of free O-H hydrogen bonds at 3500-3700  $\text{cm}^{-1}$  (sometimes missing) (smooth tongue)
  - (b) 1° Amide/Amine: 2 medium/sharp absorptions (fangs)
  - (c) 2° Amide/Amine: 1 medium/sharp absorption (missing a tooth)
  - (d) Carboxylic Acid: broad, jagged absorption (beard)
3. 3000 region - C-H stretching frequencies (sharp peaks)
  - (a) Alkene/Aromatic: C-H absorptions  $> 3000 \text{ cm}^{-1}$
  - (b) Alkane: C-H absorptions  $< 3000 \text{ cm}^{-1}$
4. Miscellaneous
  - (a) Ether: strong C–O absorption at 1000-1300  $\text{cm}^{-1}$  and lack of other absorptions
  - (b) Alkynes: C≡C absorption at 2200  $\text{cm}^{-1}$  (sharp, but can be smallish)
  - (c) Aromatic: Overtones (1650-2000  $\text{cm}^{-1}$ ) and OOP (675-900  $\text{cm}^{-1}$ ) pattern

Some functional groups (Carboxylic Acid/Amide, Alcohol/Amine, Alkenes/Aromatics and Alkanes/Ethers) can be particularly hard to differentiate between because many absorptions overlap. In those cases one must often look to secondary or even tertiary absorptions to tell the difference. Don't be afraid to nar-

## Experiment 29 IR Spectroscopy


row your results down to 1 or 2 types of molecules and then find the other spectrum that is similar and compare the two and determine which spectra is best representative of each.

### Procedure

---


1. Assign as many features as possible to each functional group.
  - (a) Indicate the order of identification
  - (b) Identification can be inclusive or exclusive (i.e. C=O present, or missing)
  - (c) Label secondary absorptions that help to verify each assignment.
2. Finally label (and circle it) your final assignment by functional group.
3. The spectra in Problems 1-12 are (a) 2,2-dimethylbutane, (b) 2-pentene (c) 2-pentyne (d) ortho-methyltoluene (e) 2-pentanol (f) butoxybutane (g) hexanal (h) 2-hexanone (i) pentanoic acid (j) 1-ethylhexyl butanoate (k) ethylamine (l) propanamide
4. Two example are provided on the next page.

## Example 1



Examining the first spectra, we first noted the carbonyl peak  $1700\text{ cm}^{-1}$  which indicates an Aldehyde, Ketone, Carboxylic Acid, Ester or Amide is present. We then looked at the  $3000\text{-}3500\text{ cm}^{-1}$  region and noted the presence of an O-H or N-H, which indicates a Carboxylic Acid or Amide. We then begin looking at the secondary peaks (in the fingerprint region) and can clearly see 3 absorptions associated with Carboxylic Acids. Lastly we note the absence of any N-H absorptions. We then conclude we have a Carboxylic Acid.

## Example 2



### *Experiment 29 IR Spectroscopy*

Examining the second spectra, we first note the lack of many features, most specifically at  $1700\text{ cm}^{-1}$  and between  $3000\text{-}3500\text{ cm}^{-1}$ . This indicates we most likely have an alkane of some sort or possibly an ether. We next looked for absorptions indicating different alkane features. The lack of intensity  $> 3000\text{ cm}^{-1}$  along with the lack of Overtones and C-H oop means we are unlikely to have a C=C therefore we don't have an alkene or aromatic. The lack of a peak near  $2100\text{-}2300\text{ cm}^{-1}$  means we don't have an alkyne either. This leaves us with the possibility of a simple alkane or ether. We not a large absorption between  $1000\text{-}1300\text{ cm}^{-1}$  which is indicative of a C-O stretch and conclude that we have an Ether.



Name: \_\_\_\_\_


Date: \_\_\_\_\_

Score: \_\_\_\_/80


## Results

Credit for all spectra in results section - Source: SDBSWeb : <http://riodb01.ibase.aist.go.jp/sdbs/> (National Institute of Advanced Industrial Science and Technology of Japan, 14 July 2008)

### Spectrum 1




### Spectrum 2




Experiment 29 IR Spectroscopy


Spectrum 3




Spectrum 4




Spectrum 5




Spectrum 6




Spectrum 7




Spectrum 8



Spectrum 9




Spectrum 10




Experiment 29 IR Spectroscopy

Spectrum 11



Spectrum 12





Name: \_\_\_\_\_

Date: \_\_\_\_\_

Score: \_\_\_\_/20

## Prelab Questions

---

1. Define the term "Functional Group Region" **AND** "Fingerprint Region".
  
2. Complete the table below by drawing a each of the molecules in the procedure section in the appropriate category below.

Name	Structure	Name	Structure
Alkane		Aldehyde	
Alkene		Ketone	
Alkyne		Carboxylic Acid	
Aromatic		Ester	
Alcohol		Amine	
Ether		Amide	

*Experiment 29 IR Spectroscopy*

Hello!